

## PICO Protein Detection Trial Kit

The **PICO Protein Detection Trial Kit** (#PICO-000090) can be stored at 4°C and -20°C (see product labels for storage conditions). Find the date of expiry on the contents of the kit.

### Further information:

- The PICO Protein Detection Trial Kit User Manual and the PICO Calculators: [www.actome.de/resources.html](http://www.actome.de/resources.html)
- Safety Data Sheets: [www.actome.de/resources.html](http://www.actome.de/resources.html)

### Notes before starting:

- Read the full PICO Protein Detection Trial Kit manual to learn more about the scientific background of a **CLC experiment** and **Relative Quantification**.
- Mix the dilutions in the 96-well plates by **gently pipetting up and down 30 times** while avoiding air bubbles.
- Perform **all centrifugation steps** with 1,000 rcf for 5 s at RT if not stated otherwise.

### Preparation of biological material:

1. Prepare the following chemicals and buffers:

<b>Additive C (5x stock)</b> add 500 µl PBS	<b>BSA (5x stock)</b> 20 mg BSA + 400 µl PBS	<b>EDTA-free Protease Inhibitor Cocktail (PIC), (25x stock)</b> dissolve 1 tablet of cOmplete Protease™ Inhibitor Cocktail (PIC) in 2 ml PBS
<b>Cell Lysis Buffer (LBT), (2x)</b> 200 µl Additive T (10x) 400 µl Additive C (5x) 80 µl PIC (25x) 200 µl Additive L 120 µl PBS	<b>Cell Lysis Buffer (LBTW), (1x)</b> 300 µl LBT (2x) 300 µl PBS	<b>Antibody Binding Control (ABC) Buffer, (1x)</b> 250 µl LBT (2x) 100 µl BSA (5x) 150 µl PBS

2. Thaw a vial of the supplied **recombinant human ErbB2/HER2 protein** and spin it down.

3. Prepare a dilution series of the protein in **LBTW** with a volume of **30 µl** each:

**Stock**  $\xrightarrow{1 \text{ in } 100}$  **S#1** (100x)  $\xrightarrow{1 \text{ in } 60}$  **S#2** (6,000x)  $\xrightarrow{1 \text{ in } 2}$  **S#3** (12,000x)  $\xrightarrow{1 \text{ in } 2}$  **S#4** (24,000x)  $\xrightarrow{1 \text{ in } 10}$  **S#5** (240,000x)

4. Calculate the volume of antibody stocks, LBT buffer, and PBS for the antibody mix (ABX) using the **Actome PICO Calculator** available online.
5. Set up the **binding reaction** in a 96-well PCR microplate (v-bottom) following the pipetting scheme (**Figure 1**).
6. Seal the plate, sonicate at full power for **1 min** and centrifuge the plate (~1,000 rcf, 30 s).
7. Incubate the plate at **4°C for 12 - 24 h**.

**Dilution and Digital PCR:**

8. Prepare the Master Mix, vortex for **10 s** and spin down.

697  $\mu$ l Ultrapure water  
284  $\mu$ l QIAcuity Probe Mix  
45.4  $\mu$ l PICO BL Probe  
45.4  $\mu$ l PICO P8 Probe  
36.3  $\mu$ l Coupling dPCR Mix

9. Calculate the required dilutions to target a **lambda of 0.15** using the **Actome PICO Calculator** available online.

10. Prepare a new 96-well plate for the dilution steps. Add the calculated amount of **PBS and 41  $\mu$ l Master Mix** to the wells shown in the dilution scheme (**Figure 2**).

11. Remove the adhesive foil of the plate containing the binding reaction and add **36  $\mu$ l PBS** to all wells (reflects the '10x pre-dilution of binding reaction' in **Figure 2**). Mix by pipetting.

12. Dilute the samples by transferring the predetermined volume (usually 1  $\mu$ l) from each sample to the corresponding wells of the dilution plate and mix by pipetting (DS 1). Repeat the dilution once more (DS 2) and finally transfer **1  $\mu$ l** to the wells containing the Master Mix (**Figure 2**).

13. Mix by pipetting and transfer **40  $\mu$ l into a QIAcuity Nanoplate 26k 24-well**. Seal the plate according to the QIAcuity user manual and run a dPCR using the following parameters:

**Priming** - QIAGEN Standard Priming Profile

**PCR conditions**


Hot-start	95°C for 2 min
Cycling	40 times
Denaturing	95°C for 15 sec
Annealing	58°C for 30 sec

**Imaging conditions**


*PICO P8 Probe* - FAM green channel, 500 ms integration time, gain 6

*PICO BL Probe* - HEX yellow channel, 400 ms integration time, gain 6

20. Instruction for evaluation of raw complex calculation with **AMULATOR** and **Relative Quantification** can be found in the full user manual.



**Figure 1:** Plate layout of the binding reaction. Five dilution steps, with four technical replicates each, were made. For the ABC three technical replicates are sufficient.



**Figure 2:** Plate setup for dPCR pre-dilution. The '10x pre-dilution of binding reaction' is prepared in the plate containing the binding reaction (see Figure 1). The two step dilution is prepared in an additional 96-well plate.

Scan the QR code for the full user manual:

