

PICO Amplification Core (AMC) Kit

The **PICO Amplification Core Kit** (#PICO-000010) contents must be stored at 4°C or -20°C for up to 6 months (see product labels for exact storage conditions). Find the date of expiry stated on the contents of the kit.

Further information:

- The PICO Handbook and the PICO AMC Kit User Manual: www.actome.de/resources.html
- Safety Data Sheets: www.actome.de/resources.html

Notes before starting:

- Perform **all centrifugation steps** with 400 rcf for 5 min at RT in a centrifuge equipped with a rotor for 15 ml falcon tubes if not stated otherwise.
- Mix the dilutions in the 96-well plates by **gently pipetting up and down 30 times** while avoiding air bubbles.
- We recommend using at least **four technical replicates** for each sample.
- We recommend using a total of **1 million cells** in a concentration of 1×10^4 cells/ μ l. However, the total input can be adapted. We also recommend using **3-fold dilution** steps for high abundance proteins and **2-fold dilution** steps for low abundance proteins.
- Sample preparation is dependent on the sample type used. Here, **cultivated cells** are described. For preparation of other samples see the full manual. Steps involving cross-linking are optional.

Start of the PICO AMC workflow (sample preparation):

1. Prepare the following chemicals and buffers:

Additive C (5x stock) add 500 μ l PBS	BSA (5x stock) 20 mg BSA + 400 μ l PBS	EDTA-free Protease Inhibitor Cocktail (PIC), (25x stock) dissolve 1 tablet of cOplete Protease™ Inhibitor Cocktail (PIC) in 2 ml PBS	PIC-PBS (1x) dissolve 1 PIC tablet in 500 μ l PBS
Cell Lysis Buffer (LBT), (2x) 200 μ l Additive T (10x) 400 μ l Additive C (5x) 80 μ l PIC (25x) 200 μ l Additive L 120 μ l PBS	Cell Lysis Buffer (LBTW), (1x) 300 μ l LBT (2x) 300 μ l PBS	Antibody Binding Control (ABC) Buffer, (1x) 250 μ l LBT (2x) 100 μ l BSA (5x) 150 μ l PBS	BS3 Cross-linker Solution (BS3S) 2 mg BS3 in 37.8 μ l DMSO BS3 Working Solution (BS3W), (1x) 6 μ l BS3S 114 μ l PBS

2. Harvest the cells with the preferred method and wash them in an appropriate volume (**3 - 15 ml**) of **PIC-PBS** by centrifugation.
3. Resuspend the cell pellet in **1 ml PBS**, transfer it to a 1.5 ml reaction tube and centrifuge again. Repeat this step once more. **IMPORTANT:** Do **not** use PIC-PBS, otherwise the PIC will quench the BS3 cross-linker.
4. Count cells in **1 ml PBS**, transfer **2×10^6 cells** into a new 1.5 ml reaction tube, centrifuge and discard supernatant.
5. Add **100 μ l of BS3W**, mix by pipetting up and down and incubate for 30 min at RT.
6. Add **900 μ l of PIC-PBS**, centrifuge and discard supernatant.
7. Count the cells again in **1 ml PBS**, transfer **1×10^6 cells** into a new 1.5 ml reaction tube, centrifuge and discard supernatant.
8. Resuspend the cells in **100 μ l LBTW** and vortex for 10 s. Incubate the cells for 3 h at 4°C.
9. Sonicate the lysate for 5 min at full power in an ultrasonic bath at room temperature.

DAY 1

10. Transfer the lysate into a **QIAshredder spin column**, centrifuge (~20,000 rcf, 2 min) and transfer the flow-through to a new 1.5 ml reaction tube.
11. Calculate the volume of antibody stocks, LBT Buffer and PBS for the antibody mix using the **Actome PICO Calculator** provided online.
12. Prepare **four technical replicates** of each sample in a 96-well PCR microplate. Combine **2 µl of sample with 2 µl of ABX** per well. Include at **least three technical replicates** of ABC (**2 µl of ABC Buffer with 2 µl of ABX**). Add **4 µl of LBTW** to an empty well as **NTC**.
13. Seal the plate with adhesive foil, sonicate for 1 min at full power in an ultrasonic bath at room temperature, centrifuge the plate (~1,000 rcf, 30 s) and incubate at 4°C overnight.

Preparing and performing of the dPCR:

14. Prepare the **Master Mix**, vortex for 10 s and spin down (~1,000 rcf, 5 s).
 - 697 µl RNase-free water
 - 284 µl QIAcuity Probe Mix
 - 45.4 µL PICO BL Probe
 - 45.4 µL PICO P8 Probe
 - 36.3 µL Coupling dPCR Mix
15. Calculate the required dilutions to target a **lambda of 0.15** using the **Actome PICO Calculator**.
16. Prepare a new 96-well plate for the dilution steps. Add the calculated amount of **PBS and 41 µl of Master Mix** to the wells shown in the dilution scheme (**Figure 1**).
17. Remove the adhesive foil of the plate containing samples and add **36 µl of PBS** to all wells. Mix by pipetting.
18. Dilute the samples by transferring **1 µl to the first set of wells** in the 96-well plate and mix by pipetting (DS1). Repeat this once more (DS 2) and finally transfer **1 µl to the wells containing the Master Mix**.
19. Mix by pipetting and transfer **40 µl into a QIAcuity Nanoplate 26k 24-well**. Seal the plate according to the QIAcuity user manual and run a dPCR using the following parameters

Priming - QIAGEN Standard Priming Profile

PCR conditions

Hot-start	95°C for 2 min
Cycling	40 times
Denaturing	95°C for 15 sec
Annealing	58°C for 30 sec

Imaging conditions

PICOact P8 Label - FAM green channel, 500 ms integration time, gain 6
PICOact BL Label - HEX yellow channel, 400 ms integrationtime, gain 6

20. Instructions for evaluation and calculators can be found in the user manual or at www.actome.de




Figure 1: Setup for dPCR pre-dilution. The '10x Pre-dilution of Binding Reaction' is prepared in a separate plate. Subsequently a two step dilution (can be adjusted if needed) is prepared in an additional 96-well plate.

Scan the QR code for the user manual:

