**BILANCIA PESAPERSONE BLUETOOTH SMART**  
**BLUETOOTH SMART BODY FAT SCALE**  
**BÁSCULA BLUETOOTH SMART**  
**PÈSE PERSONNE BLUETOOTH SMART**  
**PERSONENWAAGE BLUETOOTH SMART**  
**BALANÇA PARA PESAR PESSOAS**  
**BLUETOOTH SMART**

|                            |           |
|----------------------------|-----------|
| ISTRUZIONI PER L'USO       | <b>IT</b> |
| INSTRUCTIONS FOR USE       | <b>EN</b> |
| INSTRUCCIONES DE USO       | <b>ES</b> |
| INSTRUCTIONS D'UTILISATION | <b>FR</b> |
| GEBRAUCHSANLEITUNG         | <b>DE</b> |
| INSTRUÇÕES DE USO          | <b>PT</b> |

## BILANCIA PESAPERSONE BLUETOOTH SMART

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## BLUETOOTH SMART BODY FAT SCALE

Dear customer, **Farmamed** would like to thank you for choosing this product.

Our brand is a guarantee of high quality and highly reliable products.

For further information on this or other **Farmamed** products, please go to [www.farmamed.it](http://www.farmamed.it)

### CONTENTS

#### PAG

- 14 \_ Package contents
- 14 \_ Product description
- 15 \_ Key to symbols
- 15 \_ Safety warnings
- 15 \_ Warnings on the use of batteries
- 16 \_ General information
- 17 \_ Preparing the scale
- 17 \_ Instructions for use
- 18 \_ Assessing the results
- 21 \_ Replacing the batteries
- 21 \_ Maintenance and storage
- 21 \_ Disposal
- 22 \_ Support and warranty
- 22 \_ Troubleshooting

## PACKAGE CONTENTS

- Scale
- Instructions for use

### Functions of the appliance

This scale makes it possible to assess data on the user's fitness. The appliance is intended for private use. The scale has the following functions, which may be used by up to 8 people thanks to the appropriate memory features:

- Body weight
- Body mass index (BMI)
- Fat mass %
- Muscle mass %
- Body water %
- Visceral fat
- Skeletal system
- Basal metabolic rate


The scale also features the following functions:

- Switching of the weighing unit between kilograms "kg", pounds "lb" and stone "st"
- Automatic switch-off function
- Battery replacement indication if the batteries are low
- Automatic user recognition based on weight for using the scale offline
- Data display (7 sets of data x 2 screens)
- Weight measurement scale: 180 kg
- Division: 100 g
- LCD display

## PRODUCT DESCRIPTION


**Fig. 2: front view**

1. Display
2. Electrode pads


**Fig. 3: rear view**

3. Battery compartment
4. Unit button

**Fig.4: Display**

- |                  |                         |                             |
|------------------|-------------------------|-----------------------------|
| 1. Muscle mass % | 5. Basal metabolic rate | 9. Body weight              |
| 2. Female        | 6. Body water %         | 10. Weight readout in Kg    |
| 3. Male          | 7. Visceral fat         | 11. Weight readout in Lb    |
| 4. Fat mass %    | 8. Skeletal system      | 12. Weight readout in Stone |

## KEY TO SYMBOLS

The following symbols are used in the instructions for use:

**Warning** - Alerting to the injury or health hazard.

**Note** - Note on important information.

## SAFETY WARNINGS

Before using the appliance, carefully read the instructions for use, especially the safety warnings, and adhere to them. This manual must be considered as being part of the product and must be retained throughout its service life. Should the appliance be handed to another owner, the entire documentation should also be provided. In the event of mislaying the manual or need for more information or clarifications, contact the company at the address provided on the last page.

- Before using the product ensure the appliance is intact without any visible damage. In case of doubt do not use the appliance and contact your dealer.
- Keep the plastic bag of the packaging out of the reach of children: suffocation hazard.
- This appliance is solely intended for the use it has been designed for and in the manner indicated in the operating instructions. Any other use shall be considered misuse hence hazardous. The manufacturer shall not be liable for any damage arising from misuse or incorrect use.
- This appliance may be used independently by children no younger than 10.
- This appliance may be used by people with impaired physical, sensorial or mental abilities or by inexperienced people, only if they have been instructed on safe use beforehand and only if informed on the hazards arising from the appliance. Children must not play with the appliance. Cleaning and maintenance operations must not be performed by unsupervised children.
- In the event of fault or malfunction, switch off the appliance without tampering with it. Always contact your dealer for any repairs.
- Handle the scale with care, protect it from impact, extreme temperature changes, humidity, dust, direct sunlight and heat sources.

- DO NOT let any liquid penetrate the body of the scale. Never immerse the scale in water. Never rinse under running water. This appliance is not waterproof.
- DO NOT leave and DO NOT use the scale in water-saturated premises. Moisture exceeding 85% or contact with water or other liquids may cause corrosion and undermine proper operation of the device.
- Do not expose the appliance to impact, dust, chemical products, significant temperature changes and excessive proximity to heat sources (ovens, radiators).
- DO NOT use aggressive detergents and never immerse the appliance in water.
- DO NOT place any object on the scale when not in use.
- DO NOT step on the scale with wet feet or if its surface is wet: slipping hazard.

## WARNINGS ON THE USE OF BATTERIES

- Remove the batteries if the scale is not used for extended periods of time and store them in a cool and dry place at room temperature.
- DO NOT recharge the batteries if not rechargeable.
- DO NOT recharge the batteries (if rechargeable) with methods other than indicated in the manual provided with them or using inadequate equipment not indicated in the battery supplier documentation.
- NEVER expose the batteries to heat sources and to direct sunlight. Failure to comply with this indication may damage the batteries or cause them to explode.
- **Explosion hazard!** DO NOT throw batteries in the fire.
- Use the correct type of battery and replace them in pairs.
- The batteries must be removed or replaced by adults.
- Keep the batteries out of the reach of children: swallowing batteries is a death hazard. In the event of swallowing immediately seek medical attention.
- The acid contained in the batteries is corrosive. Avoid contact with the skin, eyes or clothes. If the battery liquid comes in contact with the skin or eyes, rinse immediately the parts concerned with water and consult a doctor.
- If the liquid leaks out of the battery, wear protective gloves and clean the battery compartment with a dry cloth. Dispose of cleaning materials.

## GENERAL INFORMATION

- The appliance is solely intended for personal use.
- Take into account that there may be measuring tolerances:  $\pm 0.3\text{kg}$  up to  $50\text{kg}$ ,  $\pm 0.4\text{kg}$  up to  $100\text{kg}$ ,  $\pm 0.5\text{kg}$  up to  $150\text{kg}$ .
- The percentage of body fat and water, skeletal muscle and the basal metabolism calculated by the scale must only be considered as indicative data. Consult a doctor for any additional information.
- Through the electrodes placed on the platform, the scale sends a weak, low intensity electrical current through the body mass.

### Recommendations for use

- In order to obtain comparable results, it is recommended to possibly weigh yourself always at the same time (preferably in the morning), before breakfast and with no clothes.
- Important during measurement: the fat mass can only be calculated with bare feet and if required, slightly moisten the soles.
- If the soles are completely dry or calloused, the results may not be satisfactory as conductivity is insufficient.
- The feet must be correctly placed on the electrodes (see picture 1).
- Keep upright and stand still during the weighing process.
- Wait a few hours before making a measurement if you have unusually exerted yourself.
- Wait about 15 minutes after getting up before making a measurement, to allow for the body water to be evenly distributed.
- Short term weight fluctuations over few days generally depend on a loss of liquids.
- Each time the appliance is moved or lifted then placed down, the first weigh-in must not be considered indicative, due to the necessary settlement of the electronic board.

Picture 1

### Limitations

It is not advisable for the following persons to use the scale because the data it measures may not be reliable:

- bearers of pace-maker or other electronic medical materials;
- very fit people or bodybuilders;
- people with fever, undergoing dialysis, with symptoms of oedema or suffering from osteoporosis;
- children under 10 years of age;
- pregnant women;
- people who take medications for the cardiovascular system;
- people who take vasodilators or vasoconstrictors;
- people with anatomical abnormalities of the legs that affect their height (significantly longer or shorter legs than average).

## PREPARING THE SCALE

### 1) INSERT BATTERIES

This scale works with three AAA 1.5V replaceable alkaline batteries (not included)

- To insert the batteries open the battery compartment by pulling the cover tab.
- Insert the batteries observing the indicated polarity. Close the cover.
- Place the scale on a flat and solid surface. A solid horizontal surface is conducive to correct measurements.
- When the battery is low a message is displayed on the screen.

### 2) SELECT THE UNIT OF MEASUREMENT

The scale is set to kg when first switched on.

If the units are other than kg when first switched on, or to select a different unit of measurement at any time (lb or st), press "UNIT" on the back of the appliance (Fig. 3), until the desired unit of measurement appears.

### 3) PLACE THE SCALE CORRECTLY

Place the scale on a flat and rigid surface. If the floor is uneven or not solid (e.g. carpet, rug, linoleum), the reliability of the measurement may be impaired.

## INSTRUCTIONS FOR USE

### 1) BODY WEIGHT MEASUREMENT

- Step on the scale with scale switched off, your feet parallel and correctly placed on the electrodes and your weight evenly distributed on the platform. Wait for the scale to compute your weight without moving. Do not lean onto anything during the measurement to prevent altering the readings.
- Step off the scale, your weight keeps being displayed for few seconds, then the scale switches off automatically.

The body composition data are displayed on the smartphone only with the "FitGo" App.

By using the scale in offline mode it is however possible to measure your weight, by following the steps below.

### 2) OPERATION WITH APP

#### List of compatible devices:

iOS from version 7.0, Android™ from version 4.3 with Bluetooth® 4.0.

iOS App

Android App


#### Installing the App and Registration

- Switch on Bluetooth® in the device settings.
- Scan the QR code or download the free "FitGo" App from App Store™ or from Google Play Store™.
- Follow the App's instructions to set up a new user with your personal information.
- Step on the platform with scale off, taking care to properly place your bare feet on the two electrodes.
- The scale displays the weight and the symbol indicating that the scale and device connection is active.
- At this stage, your device displays and stores the following data.


## 12 MEASURABLE DATA

### 1. 8 basic measurements

If the **offline weight measurement only** is carried out by a user who has previously registered on the App, and who is therefore already associated to a measurement log, the scale associates the registered user to the weight by an automatic recognition feature, therefore it is possible to view the data of the 8 basic measurements on the display even offline.

2. There are **4 additional measurements**, only available on the App, that can be activated in the settings:

- The scale switches off automatically after a few seconds if not in use.
- When using a smartphone, or a different device, before measuring the weight always ensure the Bluetooth® connection of your device is active and that the App is running and displayed on it.
- This scale can store up to 8 users in a single account and on a single device. Before making any measurement, always remember to select the user where you wish to store the data on the smartphone.
- Each user can indicate their ideal weight to build a customised plan.
- The data measurement is synchronised with the App: it is possible to monitor the data both during and after the measurement.

## ASSESSING THE RESULTS

### MEASURING PRINCIPLE

- This scale uses the BIA measurement principle, i.e. bioelectrical impedance analysis. In the space of few seconds, an imperceptible impulse, which is harmless and non hazardous, establishes the various body percentages.
- By measuring electrical resistance (impedance) and taking into account individual constants or values (age, height, sex, activity level) it is possible to establish the fat mass percentage and other body characteristics.
- Please note that the values measured by the scale approach the actual values obtained with medical tests with a deviation percentage. Only a specialist is able to perform the exact calculation of fat mass, water percentage, muscular mass and bone mass.

### BODY MASS INDEX (BMI)

The **Body Mass Index (BMI)** is a widely used parameter to obtain a general assessment of one's body weight and it is calculated by relating the person's height to their weight.

For adults aged 20 or older, the BMI is interpreted using weight classes that are the same for all ages and sexes. However, it is important to know that the body mass index does not make it possible to distinguish weight changes due to the extent of the muscular mass, bone mass, water content or fatty deposits. Therefore, one must not be deceived by the simplicity of BMI calculation: its assessment must be the responsibility of the doctor or other specialists

|             | VALUE       |
|-------------|-------------|
| UNDERWEIGHT | < 18,5      |
| NORMAL      | 18,5 - 24,9 |
| OVERWEIGHT  | 25,5 - 29,9 |
| OBESE       | > 30        |

## FAT MASS %

**Fat Mass** or FM is the total amount of lipids in the human body. It is usually expressed as a percentage of the total body mass, and consists of two components: primary fat and deposited fat. Primary fat is essential for the body's various functions and cannot therefore be done away with. The standard fat mass percentages are set out below. If the percentages shown on the scale are outside the figures shown as "normal" in the table, it is advisable to consult with your doctor.

| AGE     | WOMAN       |         |            |       | MAN         |         |            |       |
|---------|-------------|---------|------------|-------|-------------|---------|------------|-------|
|         | UNDERWEIGHT | NORMAL  | OVERWEIGHT | OBESE | UNDERWEIGHT | NORMAL  | OVERWEIGHT | OBESE |
| 10 - 39 | < 21        | 21 - 23 | 33 - 39    | > 39  | < 8         | 8 - 20  | 20 - 25    | > 25  |
| 40 - 55 | < 23        | 23 - 34 | 34 - 40    | > 40  | < 11        | 11 - 22 | 22 - 28    | > 28  |
| 56 - 80 | < 24        | 24 - 36 | 36 - 42    | > 42  | < 13        | 13 - 25 | 25 - 30    | > 30  |

## MUSCLE MASS %

There are three different types of **muscle tissue**: smooth, skeletal striated and cardiac striated.

Smooth muscle consists of smooth fibres, which are found on the walls of the digestive tract, of the respiratory tract and genitals, on the wall of venous and arterial blood vessels. In man, the skeletal striated muscle tissue forms the skeletal musculature and the muscles of the face. The cardiac striated muscle is the heart. This scale calculates the % of skeletal muscle and smooth muscle tissue in the body. Men usually have a more pronounced musculature than women and in general the percentage may increase in athletes. The muscular mass percentage generally falls within the following ranges:

|        | WOMAN | MAN   |
|--------|-------|-------|
| NORMAL | > 34% | > 40% |

## BODY WATER %

The **body water** % differs in women and men: in healthy adult women it ranges between 45% and 60%, whereas for men it is between 50% and 65%. The water percentage is higher in childhood and decreases with age. Changes in body weight may occur in the short (hours), medium (days) or long term (weeks). Small weight changes in a matter of hours are wholly normal and are due to water retention. When the weight decreases in the space of few days, it is most likely due to total water.

The water percentage generally falls within the following ranges:

| AGE     | WOMAN       |         |            | MAN         |         |            |
|---------|-------------|---------|------------|-------------|---------|------------|
|         | UNDERWEIGHT | NORMAL  | OVERWEIGHT | UNDERWEIGHT | NORMAL  | OVERWEIGHT |
| 10 - 15 | < 57        | 57 - 67 | > 67       | < 58        | 58 - 72 | > 72       |
| 16 - 30 | < 47        | 47 - 57 | > 57       | < 53        | 53 - 67 | > 67       |
| 31 - 60 | < 42        | 42 - 52 | > 52       | < 47        | 47 - 61 | > 61       |
| 61 - 80 | < 37        | 37 - 47 | > 47       | < 42        | 42 - 56 | > 56       |

Body fat does not contain much water. As a consequence, people with a high fat mass may have a lower body water percentage than standard. Whereas those who practice endurance sports, on the other hand, the percentages may be exceeded due to low fat mass and high muscle mass.

The calculation of body water with this scale does not allow medical conclusions to be drawn, e.g. as regards age-related water retention. If required consult with your doctor.

## VISCERAL FAT (V-FAT)

**Visceral or abdominal fat** indicates the accumulation of energy that is not consumed and captured in our body, which is not consumed for a variety of reasons, such as a sedentary life, lack of physical exercise, a poor diet.

The fat that is deposited at abdominal level but also in all the other areas of our body builds up due to the high presence of carbohydrates owing to excessive intake. Carbohydrates are the best source of energy for the human body but, if not properly consumed, they tend to accumulate excessively.

Abdominal fat is easily measured with a tape measure – in men normally it should not exceed 102 cm, while in women it should not exceed 88 cm.

SKELETAL SYSTEM 

The **skeleton** is the strong and flexible structure that supports our body, shaping it, protecting the internal organs, anchoring the muscles and allowing movement.

It mainly consists of bones and cartilage. Most of the bones of the female skeleton have the same shape as those of the male skeleton, albeit smaller. This scale indicates the weight of the skeletal system and this must be related to the age, height and lifestyle of each person.

BASAL METABOLIC RATE 

The **basal metabolic rate** (or BMR) is the minimum amount of calories needed by the human body to survive, therefore to perform the basic processes such as circulation, breathing, digestion, maintaining body temperature.

The basal metabolism is calculated on a person at mental and physical rest in supine position, who has been fasting for 12 hours, after at least 8 hours of sleep and with a constant ambient temperature of about 20°C. In adults, the basal metabolism ranges between 1200 and 1800 kcal. The BMR is affected by a variety of factors: age, sex, height, body composition, fever / stress, nutrition and ambient temperature.

| WEIGHT<br>Kg | WOMAN    |         |         |      | MAN  |         |         |      |
|--------------|----------|---------|---------|------|------|---------|---------|------|
|              | AGE < 17 | 18 - 23 | 30 - 50 | > 50 | < 17 | 18 - 23 | 30 - 50 | > 50 |
| < 50         | 1265     | 1180    | 1085    | 1035 | 1350 | 1200    | 1115    | 1075 |
| 50 - 65      | 1392     | 1298    | 1194    | 1139 | 1485 | 1320    | 1227    | 1183 |
| 65 - 80      | 1771     | 1652    | 1519    | 1449 | 1890 | 1680    | 1561    | 1505 |
| 80 - 95      | 2151     | 2006    | 1845    | 1760 | 2295 | 2040    | 1896    | 1828 |
| > 95         | 2530     | 2360    | 2171    | 2070 | 2700 | 2400    | 2230    | 2150 |

With weight being equal, an adolescent consumes more than an adult, and an adult more than an elderly person. What is more, always with weight being equal, men consume more than women.

OPTIONAL DATA ONLY AVAILABLE VIA THE APP

PROTEINS 

**Proteins** are an important component of human cells and their arrangement accounts for about 18% of the total mass of the human body. Muscle cells contain large amounts of protein, and it is one of the main factors that reflects a person's nutritional status, physical development and health.

DEGREE OF OBESITY 

The **degree of obesity** is the difference between the actual weight and the ideal weight.

BODY AGE 

The calculation of **body age** is tied to the basal metabolic rate. The basal metabolic rate, in turn, includes various indices including weight, body fat, muscle. Body age is a measurement of a person's physical condition, which exceeds real age. In other words, it is the criterion to assess the body's condition.

LEAN WEIGHT 

Also referred to as **lean body mass** (LBM), it is the weight of body composition after removing fat. Muscle is one of the main parts of body weight.

## TEMPORAL RELATION OF THE RESULTS

- One ought to distinguish sudden changes over a few days, medium-term variations (over weeks) and long-term variations (months).
- As a general rule, it can be considered that short-term weight variations are mainly brought about by changes in water content, while medium and long-term ones usually also involve fat and muscle mass.
- If the weight decreases in the short term while the fat mass increases or remains stable, it is probably a loss of fluids, e.g. due to training, sauna or diet only aimed at rapid weight loss.
- If the weight increases in the medium term while the fat mass decreases or remains stable, an increase in muscle mass may have occurred.
- If the weight and fat mass decrease at the same time, it means that one is losing weight.
- Fat mass, water percentage or muscle mass cannot be added (muscle tissues also contain a water percentage).

## REPLACING THE BATTERIES

When the display shows "LO" or the battery icon, replace the batteries with similar ones.

1. Open the battery compartment by pulling the cover tab.
2. Extract the spent batteries and dispose of them as indicated in the paragraph on "Disposal".
3. Insert new batteries paying attention to the indicated polarity.
4. Dispose of the batteries in accordance with local regulations.
5. Close the battery compartment.

## MAINTENANCE AND STORAGE

- Store the appliance in a cool and dry place.
- Store the appliance in a horizontal position.
- Clean the scale only with a soft and dry cloth, without adding any substances
- Never use chemical or abrasive products.

## DISPOSAL

The product packaging consists of recyclable materials. Dispose of it in compliance with environmental local regulations.

For environmental reasons, the appliance must not be disposed of with normal waste. It must be disposed of in the appropriate collection centres. Dispose of the appliance according to European Directive 2012/19/ EU on waste electrical and electronic equipment. For more information or clarifications, contact the competent local authorities for disposal.

## SUPPORT AND WARRANTY

**This appliance is solely intended for domestic use and cannot be used for professional purposes.**

- The appliance is guaranteed for 2 years from the date of purchase which must be evidenced by proof of purchase. This term complies with the laws in force ("Consumer Code" It. Legislative Decree No. 206 dated 06/09/2005) and is only applicable if the consumer is a private citizen.

- The warranty only covers manufacturing defects and is not valid if the damage is caused by an accidental event, incorrect use, negligence or misuse of the product. Excluded from the warranty are the expendable and wearing parts.

- Do not open the appliance for any reason. In the event of opening or tampering, the warranty is forfeited.

- After 2 years from purchase, the warranty expires, and technical support will be provided for a consideration. Information on technical support, be it under warranty or for a consideration, can be obtained by contacting the company at the address provided on the last page. The consumer does not owe any form of consideration for product repairs and replacements that fall under the warranty. In the event of faults contact the manufacturer.

DO NOT ship directly to Farmamed but **phone +39.0376.283211 or send an mail to [service@farmamed.it](mailto:service@farmamed.it)**. first. No warranty servicing – including to replace the product or parts thereof – shall extend the term of the original warranty of the replaced product.

- The manufacturer disclaims any liability for any damage that may arise to people, property and animals, either directly or indirectly, as a consequence of failure to comply with all the requirements set out in this documentation, especially those concerning the warnings on installation, use and maintenance of the appliance.

- Since Farmamed constantly strives to improve its products, it is entitled to modify all or part of its products without prior notice in connection with production requirements, without incurring any liability by Farmamed or its dealers.

## TROUBLESHOOTING

| Problem                                                               | Possible cause                                                                                          | Solutions                                                                                        |
|-----------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------|
| The scale does not switch on                                          | You have not stepped on the scale in the right way                                                      | Step on the scale with parallel feet. Attention: it only switches on when you are standing on it |
|                                                                       | The batteries are low                                                                                   | Replace the batteries                                                                            |
| The display shows "Lo"                                                | The batteries are low                                                                                   | Replace the batteries                                                                            |
| The display shows "Err"                                               | The weight exceeds the capacity of the scale or is less than the minimum weight measurable by the scale | Immediately remove the weight from the scale in order not to damage it and check its capacity    |
| The "FitGo" App does not show any data                                | The Bluetooth® has not been switched on in the smartphone settings                                      | Switch on Bluetooth® in the smartphone                                                           |
|                                                                       | The smartphone went into stand-by while weighing                                                        | Keep the App active on the smartphone while weighing to ensure the data are sent immediately     |
|                                                                       | You have not set up your user in the App before weighing                                                | Set up the user in the "FitGo" App                                                               |
| The "FitGo" App only displays the body weight and the body mass index | You have not stepped on the scale with bare feet                                                        | Step on the scale with bare feet                                                                 |
|                                                                       | The skin of your feet is too dry                                                                        | Slightly moisten the skin of your feet                                                           |

Other indications and solutions to other problems are available on the "FitGo" App, in the "Help" section which is accessible in the "Settings" menu.

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The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Farmamed S.r.l. is under licence.

Other trademarks and trade name are those of their respective owners.

iOS, iPhone, iPad mini, iPod touch, App Store are trademarks of Apple Inc., registered in the U.S. and other countries - Android™, Play Store™ and Google Play™ are trademarks of Google Inc.

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**Farmamed S.r.l.** – Via Cisa, 2