

UriChek[®]

2-PARAMETER

URINE REAGENT STRIPS

GLUCOSE & PROTEIN

Fast • Simple • Reliable
Urine Screening



2 ESSENTIAL PARAMETERS
Glucose & Protein



RAPID RESULTS
Results in approximately 1-2 minutes



EASY VISUAL INTERPRETATION
Compare colour development
with the provided colour chart



SIMPLE & HYGIENIC
Single-use reagent strips



CONVENIENT TESTING
Suitable for routine
urine screening



RELIABLE PERFORMANCE
Designed with minimized interference
from routine substances

TEST PARAMETERS



GLUCOSE
Helps identify glucose
in urine through an
enzymatic reaction with
visible colour development.



PROTEIN
Detects urinary protein
using the protein-error-of-
indicators principle, producing
a characteristic colour change.



IDEAL FOR



Hospitals



Diagnostic
Laboratories



Clinics



Pathology
Centres



Routine Urine
Screening



Address: F-4, 1st Floor, Plot No. 16, D-Block Tagore Nagar,
on Ajmer-Delhi, 200 Feet Bypass Rd, Jaipur,
Rajasthan 302021



Mobile:
+91 99833 33469



Email:
rajbiosis@yahoo.in



Website:
www.rajbiosis.com

QUICK URINE SCREENING

Dip. Wait. Compare.

HOW IT WORKS

- 01 **COLLECT**
Collect a fresh urine sample in a clean container.
- 02 **DIP**
Dip the reagent strip into the urine sample as instructed.
- 03 **WAIT**
Remove the strip and allow the reagent areas to develop.
- 04 **COMPARE**
Compare the colour changes with the colour chart supplied with the product.
- 05 **INTERPRET**
Determine the approximate glucose and protein level from the corresponding colour scale.

PRODUCT SPECIFICATIONS

	Product	UriChek 2-Parameter Urine Reagent Strips
	Parameters	Glucose & Protein
	Sample Type	Urine
	Test Type	Semi-quantitative
	Reading Time	Approximately 1-2 minutes
	Result Interpretation	Visual colour comparison
	Strip Type	Dry reagent strip
	Usage	Single use
	Storage	2°C – 30°C, dry place
	Humidity	Below 80% RH
	Application	Diagnostic, Clinical & Laboratory

WHY URICHEK?

- 

FAST SCREENING
Get useful urine screening information within minutes.
- 

EASY OPERATION
No complex equipment is required for visual strip testing.
- 

TWO IMPORTANT PARAMETERS
Conveniently screen for both urinary glucose and protein using one strip.
- 

PRACTICAL FOR ROUTINE TESTING
Suitable for hospitals, clinics, laboratories and routine diagnostic applications.

PRODUCT HIGHLIGHT

GLUCOSE
+
PROTEIN

Two Essential
Urine Parameters.
One Convenient
Test Strip.

Biogenix Inc. Pvt. Ltd.
www.biogenixinc.com



Address: F-4, 1st Floor, Plot No. 16, D-Block Tagore Nagar,
on Ajmer-Delhi, 200 Feet Bypass Rd, Jaipur,
Rajasthan 302021



Mobile:
+91 99833 33469



Email:
rajbiosis@yahoo.in



Website:
www.rajbiosis.com