



Raj Biosis Private Limited



FirstView™ URINE REAGENT STRIPS

Fast • Reliable • Easy-to-Use

ACCURATE RESULTS
FOR EFFECTIVE
DECISION MAKING



Quick & Convenient Testing
Results in seconds



Easy Visual Interpretation
Clear & accurate results



Multiple Parameters
Comprehensive urinalysis
in a single strip



High Accuracy & Reliability
Trusted performance for
confident diagnosis



Ideal for Routine Urinalysis
Suitable for laboratories,
clinics & hospitals



IDEAL FOR



Routine
Urinalysis



Clinical
Laboratories



Hospitals



Diagnostic
Centres



Clinics

SIMPLE TESTING.
CLEAR RESULTS.
BETTER WORKFLOW.

Raj Biosis Private Limited

F-4, 1st Floor, Plot No. 16, D-Block Tagore Nagar,
on Ajmer-Delhi, 200 Feet Bypass Rd,
Jaipur, Rajasthan 302021, India.



+91 99833 33469



rajbiosis@yahoo.in



www.rajbiosis.com

COMPREHENSIVE PARAMETERS

FirstView™ URS-11T MAU

-  Urobilinogen
-  Bilirubin
-  Ketone
-  Blood
-  Protein
-  Microalbumin
-  Nitrite
-  Leukocytes
-  Glucose
-  Specific Gravity
-  pH

FirstView™ URS-2P

-  Glucose
-  Protein



EASY TO USE – 3 SIMPLE STEPS


1



Dip

Dip the strip in fresh urine sample for 1–2 seconds.

2



Wait

Remove and shake off excess urine. Wait for the time specified.

3



Read

Compare the strip with the color chart and read the results.

PRODUCT SPECIFICATIONS

Product	FirstView™ URS-11T MAU / URS-2P
Type	Urine Reagent Strips
Storage	Store at 2–30°C. Keep away from light and moisture.
Shelf Life	24 Months
Pack Size	100 Strips / Pack
Sample Type	Urine
Method	Colorimetric
Application	For in vitro diagnostic use only



WHY CHOOSE FIRSTVIEW™?



High
Accuracy



Easy
Interpretation



Reliable
Performance



Saves Time
& Effort



Cost
Effective

IDEAL FOR



Routine
Urinalysis



Clinical
Laboratories



Hospitals



Diagnostic
Centres



Clinics

Raj Biosis Private Limited

📍 F-4, 1st Floor, Plot No. 16, D-Block Tagore Nagar,
on Ajmer-Delhi, 200 Feet Bypass Rd,
Jaipur, Rajasthan 302021, India.



+91 99833 33469



rajbiosis@yahoo.in



www.rajbiosis.com

