



Raj Biosis Private Limited

Mission®

10 PARAMETER URINALYSIS REAGENT STRIPS

FAST • CONVENIENT • COMPREHENSIVE

Mission® Urinalysis Reagent Strips provide a convenient and rapid method for qualitative and semi-quantitative analysis of multiple urine parameters. Designed for routine urinalysis, the strips help healthcare professionals screen important indicators associated with metabolic, renal, hepatic and urinary tract conditions.



10 PARAMETERS IN ONE STRIP

 1 GLUCOSE	 6 pH
 2 BILIRUBIN	 7 PROTEIN
 3 KETONE	 8 UROBILINOGEN
 4 SPECIFIC GRAVITY	 9 NITRITE
 5 BLOOD	 10 LEUKOCYTES

KEY BENEFITS

- ✓ 10 important urine parameters
- ✓ Rapid dip-and-read testing
- ✓ Easy visual interpretation using the color chart
- ✓ Qualitative & semi-quantitative results
- ✓ Convenient for routine urinalysis
- ✓ Suitable for professional laboratory and point-of-care settings
- ✓ Compatible Mission configurations are available for analyzer-based reading



IDEAL FOR

Diagnostic Laboratories • Hospitals • Clinics • Pathology Centres • Healthcare Facilities

F-4, 1st Floor, Plot No. 16, D-Block Tagore Nagar,
on Ajmer-Delhi, 200 Feet Bypass Rd,
Jaipur, Rajasthan 302021, India

+91 99833 33469

rajbiosis@yahoo.in

www.rajbiosis.com



Raj Biosis Private Limited

www.rajbiosis.com

COMPLETE URINE SCREENING MADE SIMPLE

Mission® 10 Parameter Reagent Strips combine multiple routine urine chemistry tests on a single strip, helping laboratories perform efficient screening with minimal sample handling.

TEST PARAMETERS

No.	Parameter	Helps Screen / Assess
1	GLUCOSE	Glucose abnormalities
2	BILIRUBIN	Bilirubin-related abnormalities
3	KETONE	Ketone bodies
4	SPECIFIC GRAVITY	Urine concentration
5	BLOOD	Blood / hemoglobin in urine
6	pH	Urine acidity / alkalinity
7	PROTEIN	Proteinuria screening
8	UROBILINOGEN	Urobilinogen abnormalities
9	NITRITE	Possible bacterial UTI screening
10	LEUKOCYTES	Possible urinary inflammation / UTI screening



The combination of **LEUKOCYTES** and **NITRITE** is particularly useful for screening for possible urinary tract infection, while the broader panel supports routine metabolic, kidney and liver-related screening.



SIMPLE 3-STEP TESTING

1 DIP

Immerse the reagent strip in a fresh urine specimen.



2 WAIT

Allow the reagent areas to react as per the instructions.



3 READ

Compare the colour reactions with the colour chart provided on the label.



PRODUCT HIGHLIGHTS

10

PARAMETERS



RAPID
SCREENING



EASY
HANDLING



VISUAL
READING



PROFESSIONAL
IVD USE

PACKAGING



100 TEST STRIPS
PER CANISTER

F-4, 1st Floor, Plot No. 16, D-Block Tagore Nagar,
on Ajmer-Delhi, 200 Feet Bypass Rd,
Jaipur, Rajasthan 302021, India



+91 99833 33469



rajbiosis@yahoo.in



www.rajbiosis.com