

OSCAR

URINE STRIPS 4P

Urinalysis Reagent Strips

For Semi Quantitative Detection of
Glucose, Protein, pH & Specific Gravity
in Human Urine

FAST | CONVENIENT | RELIABLE
URINE SCREENING



GLUCOSE

Screening for
glucose levels
in urine



PROTEIN

Helps detect
urinary protein



pH

Indicates the
acidity/alkalinity
of urine



SPECIFIC GRAVITY

Provides an indication
of urine concentration

Why OSCAR 4P?

- ✓ 4 essential urine parameters in one strip
- ✓ Quick and easy dipstick testing
- ✓ Convenient visual interpretation
- ✓ Suitable for routine urine screening
- ✓ Helps support efficient laboratory workflows
- ✓ Compact and easy-to-handle format

*Small Strip
Big Insights*



100

Test Strips / Pack



**For In-Vitro
Diagnostic Use Only**

Raj Biosis Private Limited

Your Partner in Better Diagnostics

F-4, 1st Floor, Plot No. 16, D-Block Tagore Nagar,
on Ajmer-Delhi, 200 Feet Bypass Rd, Jaipur, Rajasthan 302021, India.
+91 99833 33469
rajbiosis@yahoo.in
www.rajbiosis.com



Raj Biosis Private Limited

Your Partner in Better Diagnostics



VALUE YOUR LIFE

A Simple Urine Test
for a Healthier
Tomorrow

OSCAR URINE STRIPS 4P

Simple Testing. Essential Results.

OSCAR Urine Strips 4P offer a practical and reliable solution for routine urine screening by combining four commonly assessed parameters on a single reagent strip.



Reliable Results
for Routine
Urinalysis

TEST PARAMETERS

Parameter	Purpose	Clinical Relevance
 Glucose	Screening for glucose in urine	- Helps in diabetes screening - Indicates abnormal glucose levels in urine
 Protein	Helps detect urinary protein	- Useful in kidney function screening - Indicates possible renal disorders
 pH	Indicates the acidity/alkalinity of urine	- Helps assess acid-base balance - Useful in detecting urinary tract abnormalities
 Specific Gravity	Provides an indication of urine concentration	- Helps evaluate hydration status - Useful in renal function assessment

APPLICATIONS

-  Hospitals & Clinics
-  Diagnostic Laboratories
-  Pathology Centres
-  Healthcare Facilities
-  Routine Medical Screening

WIDE APPLICATIONS.
BETTER HEALTH OUTCOMES.

EASY DIPSTICK TESTING

- 1 COLLECT**

Collect a fresh urine specimen in a clean container.
- 2 DIP**

Immerse the reagent strip in the urine specimen as directed.
- 3 REMOVE**

Remove the strip and allow excess urine to drain.
- 4 READ**

Compare the reagent areas with the supplied colour chart at the specified reading time.

Always follow the manufacturer's instructions for specimen handling, testing procedure and interpretation.

PRODUCT INFORMATION

Product	OSCAR Urine Strips 4P
Type	Urinalysis Reagent Strip
Parameters	Glucose, Protein, pH, Specific Gravity
Pack Size	100 Strips
Intended Use	For In-Vitro Diagnostic Use Only



Quick



Convenient



Reliable



Trusted Worldwide

Your
Trusted Diagnostic
Partner



Better Testing
Brighter Lives

Raj Biosis Private Limited

Your Partner in Better Diagnostics

F-4, 1st Floor, Plot No. 16, D-Block Tagore Nagar,
on Ajmer-Delhi, 200 Feet Bypass Rd, Jaipur, Rajasthan 302021, India.

+91 99833 33469

rajbiosis@yahoo.in

www.rajbiosis.com