

Put infection ... where it belongs



Tonosafe™

The Tonosafe™ disposable prisms are a convenient, effective method of reducing the risk of cross infection



Convenient

Tonosafe™ avoids the need for timely and laborious disinfection without compromising the original Goldmann quality.

Safe

Developed for single use only, Tonosafe™ reliably eliminates the risk of cross infection even in the most difficult situations.

Accurate

Calibrated to Goldmann standards, Tonosafe™ delivers accurate results with all Goldmann and Perkins contact tonometers.



www.haag-streit.com

**Goldmann
Applanation Tonometry**

**HAAG-STREIT
DIAGNOSTICS**



Are you working with the Gold Standard?

Goldmann Tonometry

Goldmann Applanation Tonometry remains the Gold Standard for accurate and reproducible IOP measurement.

Make your choice

No matter if you prefer Tonosafe™ or the reusable Goldmann tonometer prism, you have made the right choice! Only these two are approved by Haag-Streit and show precision mechanics together with precision optics and therefore deliver accurate results with all Goldmann and Perkins contact tonometers.

Trustworthy

Goldmann Applanation Tonometry is proven and reproducible whether you use Tonosafe™ or the reusable Goldmann tonometer prism.

- THREE** Goldmann Tonometers
- TWO** Perfect Prisms
- ONE** Gold Standard



AT 900 D - Not for sale in the USA

HAAG-STREIT AG

Gartenstadtstrasse 10
3098 Koeniz, Switzerland
Phone +41 31 978 01 11
Fax +41 31 978 02 82
info@haag-streit.com
www.haag-streit.com



1511.7220348.02010 / 01.2013