

# 3100 Series

Modern design and ergonomic comfort make this series the economical choice.

- Contoured, 18" x 18" seat cushion
- 24" Powder-coated cast aluminum base
- Casters with polyurethane tread
- 250 lb weight capacity
- 5 year warranty

### 3125B - Operator's Stool

Dual adjustable lumbar backrest with aseptic shroud  
Pneumatic height adjustment 17.5" - 24"

### 3125BFR - Operator's Stool (high height)

Dual adjustable backrest with aseptic shroud  
Pneumatic height adjustment 20" - 26"  
Adjustable footing standard

### 3145L/R - Assistant's Stool

Adjustable body support with patented ratchet  
Left and right body support models available  
Pneumatic height adjustment 22" - 31"  
Adjustable footing standard



3125BFR

3145L/R



3125B

	seamless seat	back tilt	adjustable footing	body support	height
3125B	yes	●			17.5"-24"
3125BV	yes	●			17.5"-24"
3125BFR	yes	●	●		20"-26"
3125BFRV	yes	●	●		20"-26"
3145L	yes		●	L	22"-31"
3145LV	yes		●	L	22"-31"
3145R	yes		●	R	22"-31"
3145RV	yes		●	R	22"-31"