

Optional Accessories

External Eye Fixation Lamp
 Hard-shell Travel Case
 Wheelchair Accessible Table
 Laptop
 Computer/Monitor

General

Type of Photography	Color, Digital Red-Free, Digital Cobalt
Angle of View	45°
Magnification	2X Digital
Minimal Pupil Size	4 mm (3.3 mm SP mode)
Focus Adjustment Range	Without Compensation Lens: -10 to +15D With "-" Compensation Lens: -31 to -7D With "+" Compensation Lens: +11 to +33D
	Split-line Aligning Type Focusing
Light Source	Observation Light Source: Infrared LED Flash Light Source: White LED
Eye Fixation Lamp	Internal Eye Fixation Lamp, LED Dot Matrix, Green External Eye Fixation Lamp (Sold Separately)
Working Distance	35 mm
Sensor Resolution	20.2 Megapixels
Camera	Dedicated EOS Camera for CR-2 AF (Bundled)
Built-in Monitor	3.0 inch Color LCD Monitor
Flash Intensity	Linked to the Photography Mode Can be Set Manually
Auto Function	Auto Focus, Auto Exposure, Auto Capture, Auto Fundus
Mount Movement	Front and Back 70 mm Side to Side 100 mm Up and Down 30 mm
Chin Rest Movement	60 mm

Electrical and Environmental

Power supply	AC 100V to 240V, 50/60Hz, 1 to 0.4A
Operating environment	Temperature: 50° to 86° Fahrenheit (10° to 35° Celsius) Humidity: 30% to 90% (No Condensation) Atmospheric pressure: 700 hPa to 1060hPa

Physical Characteristics

Dimensions (H x W x D)	18.6 x 12 x 19.7 inches (473 mm x 305 mm x 500 mm)
Weight	33 pounds (15 kg)