

# Intuitive imaging, best results.

Camera Module CM 900®



## New EyeSuite Imaging Software

- Most intuitive and ergonomic operation
- Efficient image editing
- Video option always included

Digital Slit  
Lamp Imaging

 **HAAG-STREIT**  
DIAGNOSTICS

Tradition and Innovation



# The intelligent choice for optimal imaging.

- Still imaging and high resolution video
- Freeze Technology – high speed capture and data transfer
- Easy Touch – precise image control
- DFC – depth of field control
- History Trigger – never miss the perfect image


## EyeSuite Imaging:

- Intuitive and ergonomic operation
- Efficient image editing
- Appealing design

Camera Module CM 900®, exceptional imaging solution for the BD 900®, BP 900® slit lamp.


Combination of the CM 900® or EyeSuite® Imaging completes the perfect imaging system.

## Freeze Technology



Capture the image at the precise moment you press the trigger.

## History Trigger



### Rewind time and freeze the moment:

The history trigger not only captures one image but records the last few seconds in real time and lets you return to select the optimum image.

Slit Lamp BD 900®



Slit Lamp BP 900®



From the leader  
in slit lamp imaging

**HAAG-STREIT AG**  
Gartenstadtstrasse 10  
3098 Koeniz, Switzerland  
Phone +41 31 978 01 11  
Fax +41 31 978 02 82  
info@haag-streit.com  
www.haag-streit.com

 **HAAG-STREIT**  
DIAGNOSTICS

Tradition and Innovation