



# CR-2 PLUS AF

Digital Retinal Camera

*The new CR-2 PLUS AF Digital Retinal Camera provides 10 important technological enhancements to the CR-2 PLUS model.*

These features include:

**1. Auto-Fundus**

Allows automatic switching from the external eye to retinal observation mode

**2. Auto-Focus**

Automatically focuses the eye when the joystick is partially pressed

**3. Auto-Capture**

Automatically captures the image when the eye is properly focused

**4. Image Error Detection**

Intelligent software automatically confirms correct alignment and focus

**5. Enhanced Visibility During Alignment**

Two times bigger/three times brighter working dots for easy patient alignment

**6. Manual Alignment Override**

Easily switch to manual focusing on demand

**7. Quick Preview**

Image previewed with EOS camera technology at the time of capture

**8. Anterior Segment Compensator**

Integrated diopter compensation on focus ring simplifies anterior photos

**9. 3-Color Image Profiles**

Selector allows clinician to determine preferred color palette

**10. External Monitor Port**

Added port accommodates external monitor via HDMI Cable



## Benefits

- Novice-friendly design
- Improved exam throughput\*
- Reduced training time\*
- Enhanced ease of use\*
- Increased image acquisition quality\*
- Capture software seamlessly integrates into imageSPECTRUM Image Management Software (sold separately)

\* As compared to the CR-2 PLUS Digital Retinal Camera



1-800-970-7227  
eccsales@usa.canon.com  
www.usa.canon.com

Healthcare Solutions Division  
Business Imaging Solutions Group  
Canon U.S.A., Inc.  
15955 Alton Parkway  
Irvine, CA 92618

CANON and EOS are registered trademarks of Canon Inc. in the United States and may also be registered trademarks or trademarks in other countries. All other referenced product names and marks are trademarks of their respective owners. Specifications and availability subject to change without notice.  
©2013 Canon U.S.A., Inc. All rights reserved.