

MARCO *idoc*® SYSTEM REQUIREMENTS

HARDWARE & SOFTWARE SPECIFICATIONS:

- OS** Windows XP 32bit / 64bit (All editions) SP2/SP3
Windows Vista 32bit / 64bit (All editions)
Windows 7 32bit / 64bit (All editions) SP1
- CPU** Intel 1 (GHz) or faster 32-bit (x86) or 64-bit (x64) processor
- RAM** 1 (GB) RAM (32-bit) or 2 GB RAM (64-bit) or more
- HD** 200 (GB) or more
- Drive** CD-RW
- (NEW ITEM) Video Card** Any video card which has enough capacity to show high quality images. (Consult your IT professional.)
- I/F Desktop** IEEE 1394A (Firewire) PCI Board is required to be installed, 1 RS232-C Serial Port (or Serial to USB adapter), 4 or more open USB 2.0 ports at the rear of the computer case.

I/F Laptop Externally powered IEEE 1394A (Firewire) card for the specific Laptop auxiliary card slot.

CAMERA SPECIFICATIONS:

- Image Device** ½" Type progressive scan SONY IT CCD
- Effective Picture Elements** Up to 1388 x 1038 pixel
- Transfer Rate** 100/200/400/800 Mb/s
- Frame Rate** 15MHz at full resolution (1388 x 1038)
- Resolution Depth** 8 bit
- Color Modes** Raw8 (Mono8), YUV4:2:2, YUV4:1:1
- Digital Interface** IEEE1394; DCAM V1.30
- Power Requirements** DC 8V-36V via IEEE1394 Cable
- Power Consumption** Less than 3 Watt(@12 VDC)

