

# Go for the Gold Standard

Applanation Tonometer AT 900®



Goldmann  
Applanation Tonometry

**HS** HAAG-STREIT  
INTERNATIONAL

Tradition and Innovation

# The Goldmann applanation tonometer


- Gold standard in IOP measurement
- Best reproducibility
- In combination with any slit lamp
- Reliable and durable
- Reusable tonometer prisms
- Tonosafe™ disposable prisms

Goldmann applanation tonometry remains the gold standard for fast and reproducible IOP measurement.


It is based on the Goldmann principle – a principle proven over more than 50 years and with millions of patients.

The positive experience of many thousands of tonometer users highlights the trustworthiness of your measurement results.


Goldmann applanation tonometry is economic, simple and reproducible. The HAAG-STREIT high precision manufacture guarantees continuous, reliable performance and ensures many years of usage without the necessity of manufacturer maintenance.




Model AT 870  
for Zeiss-type slit lamps



AT 900® model R



AT 900® model T



AT 900® model BQ

Tonometer model			
Slit lamp model	Mod. R	Mod. T	Mod. BQ
BD 900® / BC 900® *	●		
BM 900® / BMV 900® *	●	●	
BP 900®		●	●
BQ 900®		●	●

\* Discontinued slit lamp model

**From the original Goldmann manufacturer**

## HAAG-STREIT AG

Gartenstadtstrasse 10  
3098 Koeniz, Switzerland  
Phone +41 31 978 01 11  
Fax +41 31 978 02 82  
info@haag-streit.com  
www.haag-streit.com

**HAAG-STREIT**  
INTERNATIONAL

Tradition and Innovation