

Easy use, optimal imaging

Imaging Module IM 900®



New EyeSuite Imaging Software

- Most intuitive and ergonomic operation
- Efficient image editing
- Video option always included

Digital Slit
Lamp Imaging

 **HAAG-STREIT**
DIAGNOSTICS

Tradition and Innovation



The intelligent choice for optimal imaging.

- Still imaging and high resolution video
- Freeze Technology – high speed capture and data transfer
- Easy Touch – precise image control
- DFC – depth of field control
- History Trigger – never miss the perfect image

EyeSuite Imaging:


- Intuitive and ergonomic operation
- Efficient image editing
- Appealing design

The Imaging Module IM 900® is the fully integrated compact imaging solution for the BQ 900® slit lamp.

The Combination of the IM 900® with EyeSuite™ Imaging completes the perfect imaging system.




Freeze Technology



Capture the image at the precise moment you press the trigger.

History Trigger



Rewind time and freeze the moment:

The history trigger not only captures one image but records the last few seconds in real time and lets you return to select the optimum image.

Easy Touch



Control the image:

Easy touch more than just a trigger, allowing simple and ergonomic management of camera settings.

Depth of field control (DFC)



True aperture:

An unique feature of the Haag-Streit imaging system, facilitating complete control of the depth of field.

From the leader in slit lamp imaging

HAAG-STREIT AG

Gartenstadtstrasse 10
3098 Koeniz, Switzerland
Phone +41 31 978 01 11
Fax +41 31 978 02 82
info@haag-streit.com
www.haag-streit.com



Tradition and Innovation