

Welch Allyn Spot™ Binocular Autorefractor

Eye care professionals are challenged with maintaining practice flow when treating patients with limited mobility, disabilities, and children of all ages. Desktop autorefractors are confined to a single location, which can limit patient access while dictating a practice flow that may not be consistent with today's eye care practices' needs. These challenges of supporting accessible and affordable care for patients of all ages can challenge a staff's ability to maintain a high level of customer service in a highly competitive industry.



Spot™ is a new generation of handheld binocular autorefraction. A fast and unobtrusive approach to autorefraction brings versatility to your practice. Whether practice needs dictate supporting ADA requirements, capturing ACA patients, or to simply streamline practice flow, the Spot Autorefractor is an accurate, cost-effective exam support tool for today's eye care professional.

Expand Your Capabilities

- From infants to seniors, see patients as young as 6 months
- Ideal for patients with limited mobility and other disabilities
- Requires minimal patient cooperation
- As simple as taking a photo

Accurately Measures for:

- Refractions – monocular or binocular (SE, DS, DC, axis, +/- cyl conventions)
- Measurement range: -7.50D to +7.50D spherical equivalent
- Eye alignment (measure vertical and horizontal in degrees)
- Pupillary distance
- Pupil size

Easy to Implement and Use

- Quick and noninvasive
- Results of both eyes displayed in seconds
- Goes where you go – use anywhere, in the office or at events
- WiFi connectivity for easy data transfer and printing

Consistent, Objective Results

- Results displayed on screen in seconds
- Wireless printing direct from device
- Wireless label printing

Spot™ Binocular Autorefractor

Since the Spot Binocular Autorefractor is portable, you can bring it to patients with limited mobility. Therefore, in the United States, use of Spot may qualify for a 50% tax credit under the Americans with Disabilities Act (ADA). Other federal and state tax incentives could apply. These can add up to a substantial savings on your investment in a Spot Binocular Autorefractor. Consult with your tax advisor to determine if use of Spot in your practice is eligible for these tax incentives.



“ I love Spot. I take it between my two offices and use it on many patients per day. It's very fast and unobtrusive ...it's priceless. ”

—Dr. Bill Kegerize, Optometrist
Genoa Family Eye Care, Oregon, OH

Ordering Information

- VS100-B Spot Binocular Autorefractor, power supply and cord, wrist strap, and DFU
- VS100S-B Spot Binocular Autorefractor, power supply & cord, wrist strap, DFU, and carry case
- 106144 Spot Binocular Autorefractor Carry Case

To learn how Spot can help you test all of your patients efficiently and accurately, contact your local distributor, call 1.800.535.6663 to speak with a Welch Allyn customer care representative, or visit welchallyn.com/spotvision

Welch Allyn Corporate Headquarters
4341 State Street Road, P.O. Box 220
Skaneateles Falls, NY 13153-0220 USA
(p) 800.535.6663 (f) 315.685.3361

WWW.WELCHALLYN.COM

WelchAllyn
Advancing Frontline Care™