

TF-12 TRIAL FRAME

PD Range: 48-80mm



CE FDA



TF-10 TRIAL FRAME

PD Range: 48-80mm



CE FDA

TF-19 TRIAL FRAME

PD Range: 54-70mm

TITANIUM



CE FDA

TF-32 TRIAL FRAME

FIXED PD



CE FDA

PD range(10 options): 52mm, 54mm, 56mm, 58mm, 60mm;
62mm, 64mm, 66mm, 68mm, 70mm
Color available: Blue, Pink, Grey



TF-13 TRIAL FRAME FOR KIDS

FIXED PD

PD range: 48mm and 56mm are available



CE FDA

TF-23 TRIAL FRAME

FIXED PD

PD range: 56mm, 58mm, 60mm, 62mm and 64mm are available
Net weight: 40g



CE FDA

TC-1 TRIAL CLIP



CE FDA



It has spring piece, hold the trial lens in.

CP-9B PH METER

Measurer of Pupillary Height

FEATURES

Convenience, visualization, easy operation and accuracy make it indispensable during the Optometry.

INSTRUCTIONS

- 1) Clip this measurer on the frame, aim the center point at the intermediate position of the bridge.
- 2) Adjust the measuring pieces by the knobs on both sides. When the upper edge of the measuring pieces aligned with the pupil center, the pupil height can be read at the lower edge of the frame.

CP-9 PD & PH METER

Measurer of Pupillary Height and Pupillary Distance

FEATURES

Including all the functions of CP-9B, meanwhile CP-9 can measure the pupillary distance.

